

F-OCTANE LIGHT



More Information:
www.fluoron.de

Biocompatible

**Physiologically and
chemically inert**

**Significantly reduced
density of 1.46 g/cm³**

F-OCTANE^{LIGHT} is a highly pure mixture of Perfluorooctane and F6H6[®]. Its specially adjusted density promotes "gentle" reattachment of the retina.

- Intraoperative adjuvant for retinal attachment.
- The reduced density of 1.46 g/cm³ allows tractions to remain visible to the surgeon.

G-80705 F-OCTANE^{LIGHT}
Vial 5 ml, **sterile**

G-80707 F-OCTANE^{LIGHT}
Vial 7 ml, **sterile**