

# Refraction unit Smile



Refraction unit

**Med**  
**instrus**


# Refraction unit SMILE

## Optional

- Phoroptor arm
- Chin rest for slit lamp / ophthalmometer
- Foot rest for the chair
- Foot switch control for up / down movement of the chair
- Left handed version
- ~110V power supply

## Refraction unit SMILE

- The vertical movement of the chair is adjustable from integrated touch panel on the tabletop and / or from the foot-switch (optional)
- Tabletop for two instruments
- Column with halogen overhead lamp with adjustable illumination and support for chart projector
- Integrated adjustable power transformer for slit lamp / ophthalmometer
- 230V socket for additional instrument (e.g. auto lensmeter, auto phoroptor, etc.)



## Technical information

Voltage	~230V 50 / 60Hz
Total weight	160 kg
Max height	181 cm
Table height from the ground	86 cm
Halogen lamp	20W 230V
Adjustable voltage supply for slit lamp / ophthalmometer	6V and 12V
Chair vertical movement	15 cm
Chair rotation	360°
Chair height from the ground	52 cm



**UAB "Medinstrus"**  
Respublikos 20, LT - 89224 Mažeikiai,  
Lietuva (Lithuania).

**Tel:** +370 443 25451  
**Fax:** +370 443 25455  
**E-mail:** info@medinstrus.lt  
**Web:** www.medinstrus.lt

